# Pen type, all in one, IP-67

# **CONDUCTIVITY METER**

Model: PCD-431



Your purchase of this CONDUCTIVITY **MFTFR** marks a step forward for you into the field precision measurement. Although this METER is a complex and delicate instrument. its durable structure will allow many years of use if proper operating techniques are developed. Please read the following instructions carefully and always keep this manual within easy reach.



# **OPERATION MANUAL**

# **TABLE OF CONTENTS**

1. FEATURES	.1
2. SPECIFICATIONS	.2
3. FRONT PANEL DESCRIPTIONS	.5 .5 .5 .5 .5
4. MEASURING PROCEDURE	6 7 8 8
5. CALIBRATION	9
6. REPLACEMENT OF BATTERY	10
7. OPTIONAL ACCESSORIES	10

# 1. FEATURES

- \* All in one pen type conductivity meter provides fast, accurate readings with digital reading.
- \* Conductivity measurement ( uS, mS ) or TDS ( Total Dissolved Solids, ppm ) can be selected.
- \* Conductivity: Two ranges, 2,000 uS, 20.00 mS.
- \* TDS: Two ranges, 20,000 ppm, 2,000 ppm.
- \* Carbon rod electrode for long life.
- \* Build in temperature sensor, ATC ( auto temperature compensation ).
- \* Temperature measurement, °C, °F.
- \* IP67, water proof and protection.
- \* LCD with two displays show conductivity and Temp. value at same time.
- \* Data hold function for freezing the desired value.
- \* Auto power off to save the battery life.
- \* Records max. and min. value with recall.
- \* Microcomputer circuit, intelligent function, high accuracy.
- \* Compact size, light weight.
- \* Power supply by DC 1.5 V battery ( UM4/AAA ) x 4 PCs,
- \* Available for wide applications, such as aquarium, beverage, fish hatcheries, food processing, photography, laboratory, quality control, school & colleges, swimming pools, water conditions.

# 2. SPECIFICATIONS

Display	LCD, size: 20 mm x 28 mm.		
Measurement	* Conductivity ( uS, mS )		
	* TDS ( Total Dissolved Solids, ppm )		
	* Temperature (°C,°F)		
Ranges	Conductivity:		
* two ranges	2000 uS, 20.00 mS		
* auto range	TDS:		
	2,000 ppm 20,000 ppm		
Accuracy	± (2% FS + 1 d)		
* 23 ± 5 °C	* FS : full scale		
Temperature	Automatic from 0 to 60 $^{\circ}$ C (32 - 140 $^{\circ}$ F),		
Compensation	with temperature compensation factor		
	variable between 0 to 5.0% per C.		
Conductivity	Carbon rod electrode for long life.		
Probe			
Structure			
Data Hold	Freeze the display reading.		
Memory Recall	Maximum & Minimum value.		
Sampling	Approx. 0.8 second.		
Time			
Circuit	Custom one-chip of microprocessor LSI		
	circuit.		
Power off	Auto shut off saves battery life or		
	manual off by push button.		
	* Power will off automatically after		
	10 min., if no button be pressed.		
Operating	0 to 50 $^{\circ}\mathrm{C}$ ( 32 to 140 $^{\circ}\mathrm{F}$ ).		
Temperature			
Operating	Less than 80% RH.		
Humidity			

Power Supply	DC 1.5V battery ( UM-4/AAA ) x 4 PCs.	
Power	Approx. 5.7 mA.	
Consumption		
Dimension	190 x 40 x 40 mm	
	(7.5 x 1.6 x 1.6 inch).	
Weight	171 g/0.38 LB.	
Standard	Instruction Manual 1 PC	
Accessories		
Power off	Auto shut off saves battery life or	
	manual off by push button.	
	* Power will off automatically after	
	10 min., if no button be pressed.	
Standard	Instruction Manual 1 PC	
Accessories		
Optional	* Soft carrying case with sash	
Accessories	( 210 x 80 x 50 mm ), Model : CA-52A	
	* Hard carrying case	
	( 280 x 195 x 65 mm ), Model : CA-06	
	* 1.413 mS Conductivity Standard	
	Solution, Model: CD-14A	

### 2-2 Electrical Specifications (23 ± 5 $^{\circ}$ C)

A. Conductivity

Range	Measurement	Resolution	Accuracy
2000 uS	0 to 2000 uS	1 uS	± (3% F.S.+1d)
20 mS	2 to 20.00 mS	0.01 mS	* F.S Full scale

<sup>\*</sup> Auto range.

compensation factor variable between 0 to 5.0% per C.

B. TDS (Total Dissolved Solids)

Range	Measurement	Resolution	Accuracy
2,000 ppm	0 to 1,320 ppm	1 ppm	± (3% F.S.+1d)
20,000 ppm	1,320 to 13,200 ppm	10 ppm	* F.S Full scale

<sup>\*</sup> Auto range.

Automatic from 0 to 60  $^\circ$ C ( 32 - 140  $^\circ$ F ), with temperature compensation factor variable between 0 to 5.0% per  $^\circ$ C.

### C. Temperature

Function	Measuring Range	Resolution	Accuracy
$^{\circ}$	0 ℃ to 60 ℃	0.1 ℃	0.8 ℃
°F	32 °F to 140 °F	0.1 °F	1.5 °F
* @ 23± 5℃			

<sup>\*</sup> Above specification tests under the environment RF Field Strength less than 3 V/M & frequency less than 30 MHz only.


<sup>\*</sup> Temperature Compensation : Automatic from 0 to 60  $^\circ$  ( 32 - 140  $^\circ$ F ), with temperature

<sup>\*</sup> mS - milli Simens \* uS - micro Simens \*  $23\pm5\%$ 

<sup>\*</sup> Temperature Compensation :

<sup>\*</sup> ppm - parts per million \* 23 $\pm$  5 $^{\circ}$ C

# 3. FRONT PANEL DESCRIPTION



- 3-1 Battery compartment/Cover
- 3-2 Display
- 3-3 Power Button
- 3-4 Hold Button, Unit Button
- 3-5 REC Button (  $^{\circ}$ C/ $^{\circ}$ F Button )
- 3-6 Conductivity Electrode
- 3-7 Protection Cover

# 4. MEASURING PROCEDURE

#### 4-1 Conductivity measurement


- 1) Turn on the meter by pressing the "Power Button" (3-3, Fig. 1) momentarily.
  - \* Press the "Power Button" (3-3, Fig. 1) momentarily again will turn off the meter.

The "Display" (3-2, Fig. 1) will show the unit as:



Now the meter is ready for the conductivity measurement.

- 2) Hold the meter by hand and immerse the "Sensing Electrode" (3-6, Fig. 1) immersed wholly into the measured solution. Shake the "Sensing Electrode" to let the electrode's internal air bubble drift out from the sensing Electrode.
  - " Display " (3-2, Fig. 1) will show the conductivity mS (uS) values, at the same time the left bottom display will show the Temp. value of the measured solution.



#### 4-2 TDS measurement


Same procedures as above " 4-1 Conductivity measurement " but should select the measuring unit to " ppm " at first, the procedures are following:

Press " Unit Button " ( 3-4, Fig. 1 ) continuously at least two seconds until the " Display " show the unit " ppm " instead of " uS ", then release the button.

Now the meter is ready for TDS ( Total Dissolved Solids ) measurement.

#### Remark:

If the measurement value > 9990 ppm, then the last digit will present in the low display as :



#### 4-3 Data Hold

- During the measurement, press the "Hold Button" (3-4, Fig. 1) momentarily to hold the measured value. The LCD will show a "HOLD" symbol.
- 2) Press the " Hold Button " once again to release the data hold function.

### 4-4 Data Record (Max., Min. reading)

- The data record function records the maximum and minimum readings. Press the "REC Button" (3-5, Fig. 1) momentarily to start the Data Record function, shows "REC" on the display.
- 2) With the "REC" symbol on the display.
  - a) Press the "REC Button" (3-5, Fig. 1) momentarily, the "REC MAX" symbol along with the maximum value will appear on the display.
  - b) Press the "REC Button" (3-5, Fig. 1) momentarily again, the "REC MIN" symbol along with the minimum value will appear on the display.
    - \* When display shows " REC MAX " or " REC MIN ", press the " Hold Button " ( 3-4, Fig. 1 ) momentarily will delete the max. ( min. ) value, the display will show the " REC. " only and execute the memory function continuously.
  - c) To exit the memory record function, press the " REC " button for 2 seconds at least. The display will revert to the current reading, not show " REC " indicator.

# 4-5 Change Temp. unit $\mathcal{C}$ , $\mathcal{F}$

- The temperature units can be selected to °C or °F.
- 2) Pressing the "REC Button,  $^{\circ}$ C/ $^{\circ}$ F Button " ( 3-5, Fig. 1 ) continuously at least 2 seconds, then Temp. unit will change from  $^{\circ}$ C to  $^{\circ}$ F or  $^{\circ}$ F to  $^{\circ}$ C.
- After the temperature unit is selected, it will be memorized into the memory circuit. If power off and on again will present the existing selection Temp. unit.

#### 4-6 Auto power off management

The meter has built-in "Auto Power Shut-off" function in order to prolong battery life. The meter will switch off automatically if none of the buttons are pressed within 10 min..

# 5. CALIBRATION

The meter already made the calibration precisely during the production in the factory. When the user get the meter, the calibration are not necessary. However after use the meter for a long period, if necessary to make the calibration, please contact the distributor and should calibrated by the authorized technician people only.

# 6. REPLACEMENT OF BATTERY

- 1) When the LCD display show " \(\subset \subset \subs
- 2) To replace the battery, rotate and remove the "Battery Cover" (3-1, Fig. 1), take out the old baDC 1.5V battery (UM-4/AAA) x 4 PCs.
- 3) When install the batteries, should make attention the battery polarity.
- 4) After install the batteries, reinstall the battery cover again.

# 7. OPTIONAL ACCESSORIES

Hard carrying case.	Soft carrying case.
Model: CA-06	Model: CA-52A
280 x195 x 65 mm	210 x 80 x 50 mm